



Product Datasheet

CRYAA antibody (orb178603)


Description	Rabbit polyclonal antibody to CRYAA.
Reactivity	Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human CRYAA(NP_000385.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	19KDa
Uniprot ID	P02489
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of mouse kidney, liver, and heart.



Immunohistochemical staining of rat kidney tissue.



Immunohistochemical staining of human cornea tissue.