


Product Datasheet

Biotin Antibody (orb1786031)

Description	Biotin Antibody
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IP, MC
Immunogen	biotinylated limulus hemocyanin
Target	Biotin
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Isotype	Mouse IgG2a kappa
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.



ELISA
analysis of a
biotinylated
antigen...