


Product Datasheet

Human Protogenin ELISA Kit (orb1786034)

Description	Human Protogenin ELISA Kit (96 Tests). Quantitate Human Protogenin in cell culture supernatants, ser...
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: F36-A952
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA)
Uniprot ID	Q2VWP7
Sensitivity	8 pg/ml
Expiration Date	Please enquire.



Human
Protogenin
PicoKine
ELISA Kit
stan...