



Product Datasheet

ATP1B1 antibody (orb178604)


Description	Rabbit polyclonal antibody to ATP1B1
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human ATP1B1(NP_001668.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	55kDa
Uniprot ID	P05026
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human stomach tissue.



Immunohistochemical staining of human gastric tissue.



Immunohistochemical staining of mouse lung tissue.