

Product Datasheet

Mouse LIGHT / TNFSF14 Protein, His Tag, active trimer, low endotoxin (MALS verified) (orb1786126)

Description

Mouse LIGHT Protein, His Tag is

(HEK293). It contains AA Asp 72 - Val...

Reactivity

Mouse

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

LIGHT

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

N-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 58.9 kDa. The protein migrates as 55-60 kDa and 65-80 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Mouse LIGHT Protein, His Tag (orb1786126) is expressed from human 293 cells (HEK293). It contains AA Asp 72 - Val 239 (Accession # Q9QYH9-1).

MW

58.9 kDa

Uniprot ID

Q9QYH9

Expiration Date

6 months from date of receipt.