

Product Datasheet

Human CCL8 Protein, His Tag (orb1786138)

Description	Human CCL8 Protein, His Tag is expressed from human 293 cells (HEK293). It contains AA Gln 24 - Pro ...
Reactivity	Human
Conjugation	Unconjugated
Target	CCL8
Preservatives	Powder
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	N-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 10.8 kDa.
Source	Human CCL8 Protein, His Tag (orb1786138) is expressed from human 293 cells (HEK293). It contains AA Gln 24 - Pro 99 (Accession # P80075).
MW	10.8 kDa
Uniprot ID	P80075
Expiration Date	6 months from date of receipt.