

Product Datasheet

Human BMP-5 Protein, Tag Free [Out of stock] (orb1786140)

Description	This gene encodes a secreted ligand of the TGF-beta (transforming growth factor-beta) superfamily of...
Conjugation	Unconjugated
Target	BMP-5
Form/Appearance	Powder
Storage	● -20°C to -70°C for 12 months in lyophilized state; ● -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 16.6 kDa.
Source	Human BMP-5 Protein, Tag Free (orb1786140) is expressed from CHO cells. It contains AA Ala 317 - His 454 (Accession # P22003-1).
MW	16.6 kDa
Uniprot ID	P22003
Expiration Date	12 months from date of receipt.