

Product Datasheet

Human Laminin 521 Protein, His Tag [Out of stock] (orb1786146)

Description	Laminins are a major component of basement membranes. Each laminin molecule is a heterotrimeric glyco...
Conjugation	Unconjugated
Target	Laminin 521
Form/Appearance	Powder
Storage	● -20°C to -70°C for 12 months in lyophilized state; ● -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	N-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 400 kDa 194.3 kDa 176 kDa.
Source	Human Laminin 521 Protein, His Tag (orb1786146) is expressed from human 293 cells (HEK293). It contains AA Arg 36 - Ala 3695 & Gln 33 - Gln 1798 & Gln 34 - Pro 1609 (Accession # O15230-1 & P55268-1 & P11047-1).
MW	400 kDa 194.3 kDa 176 kDa
Uniprot ID	P11047, O15230, P55268
Expiration Date	12 months from date of receipt.