

Product Datasheet

**Mouse IgE Fc Protein, His Tag (MALS & SPR
verified) (orb1786148)**

Description

As one of the five designated

immunoglobulin E (IgE) plays a major role in ...

Reactivity

Mouse

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

IgE Fc

Preservatives

PBS, pH 7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 39.3 kDa. The protein migrates as 45-55 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Mouse IgE Fc Protein, His Tag (orb1786148) is expressed from human 293 cells (HEK293). It contains AA Val 91 - Ser 421 (Accession # P06336).

MW

39.3 kDa

Uniprot ID

P06336

Expiration Date

6 months from date of receipt.