

Product Datasheet

**Human Apolipoprotein E / APOE3 Protein, Tag Free
(MALS verified) (orb1786200)**

Description

APOE is an apolipoprotein, a

particles, that mainly functions in lipo...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Apolipoprotein E

Preservatives

PBS, pH 7.4

Form/Appearance

Powder

Storage

-20°C

Tag

Native

Note

For research use only

Application notes

This protein carries no "tag". The protein has a calculated MW of 35.0 kDa. The protein migrates as 37-40 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Human APOE3 Protein, Tag Free (orb1786200) is expressed from human 293 cells (HEK293). It contains AA Lys 19 - His 317 (Accession # P02649-1).

MW

35.0 kDa

Expiration Date

6 months from date of receipt.