


Product Datasheet


ANP32A antibody (orb178621)

Description	Rabbit polyclonal antibody to ANP32A
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human ANP32A(NP_006296.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	28kDa
Uniprot ID	P39687
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.




Immunohistochemical staining of human B

...



Immunohistochemical staining of Rat liver...



Immunohistochemical staining of human gastric...