



Product Datasheet

CXXC1 antibody (orb178657)


Description	Rabbit polyclonal antibody to CXXC1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human CXXC1(NP_001095124.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	75kDa
Uniprot ID	Q9P0U4
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200, ChIP: 1:20-100
Expiration Date	12 months from date of receipt.



WB analysis of extracts of various cell ...



Immunofluorescence analysis of A549 cell...



Immunohistochemical staining of mouse lu...