


Product Datasheet

Human Flt3 Ligand ELISA Kit (orb178672)

Description	Human Flt3 Ligand ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: T27-P185
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	26416 MW
Uniprot ID	P49771
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human Flt-3
ligand
PicoKine
ELISA Kit
st...