


Product Datasheet

Mouse Flt-3 Ligand ELISA Kit (orb178673)

Description	Mouse Flt-3 Ligand ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: G27-R188
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
MW	26141 MW
Uniprot ID	P49772
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:5, concentration in serum is around 1000 pg/ml.
Expiration Date	Please enquire.



Mouse Flt-3 Ligand PicoKine ELISA Kit sta...