


Product Datasheet

Rat Fractalkine/CX3CL1 ELISA Kit (orb178674)

Description	Rat Fractalkine/CX3CL1 ELISA Kit
Reactivity	Rat
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q25-T334
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
MW	41977 MW
Uniprot ID	055145
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Rat
Fractalkine
PicoKine
ELISA Kit
stand...