


Product Datasheet

Human IGFBP2 ELISA Kit (orb178675)

Description	Human IGFBP2 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E40-Q328
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and urine.
MW	34814 MW
Uniprot ID	P18065
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:200, concentration in serum and plasma is around 200 ng/ml.
Expiration Date	Please enquire.



Human
IGFBP2
PicoKine
ELISA Kit
standard...