


Product Datasheet

Mouse VEGFR2/KDR ELISA Kit (orb178680)

Description	Mouse VEGFR2/KDR ELISA Kit
Reactivity	Mouse
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A20-E762
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
MW	152517 MW
Uniprot ID	P35918
Sensitivity	2 pg/ml
Dilution Range	Dilution ratio of 1:4, concentration in serum is 20-40 ng/ml.
Expiration Date	Please enquire.



Mouse
VEGFR2
PicoKine
ELISA Kit
standard...