


## Product Datasheet

### Human CD36/SR-B3/GP4 ELISA Kit (orb178683)

|                            |                                                                                                                                                                                   |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Human CD36/SR-B3/GP4 ELISA Kit                                                                                                                                                    |
| <b>Reactivity</b>          | Human                                                                                                                                                                             |
| <b>Range</b>               | 156 pg/ml - 10,000 pg/ml                                                                                                                                                          |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                             |
| <b>Immunogen</b>           | Expression system for standard: NS0; Immunogen sequence: G30-N439                                                                                                                 |
| <b>Storage</b>             | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| <b>Note</b>                | For research use only                                                                                                                                                             |
| <b>Sample Types</b>        | cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).                                                                                                        |
| <b>MW</b>                  | 53053 MW                                                                                                                                                                          |
| <b>Uniprot ID</b>          | <b>P16671</b>                                                                                                                                                                     |
| <b>Sensitivity</b>         | 10 pg/ml                                                                                                                                                                          |
| <b>Dilution Range</b>      | Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 156 pg/ml.                                                                             |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                   |



Human  
CD36/SR-B3  
PicoKine  
ELISA Kit  
stan...