


Product Datasheet

Human CXCL7/PPBP/NAP-2 ELISA Kit (orb178686)

Description	Human CXCL7/PPBP/NAP-2 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A59-D128
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	13894 MW
Uniprot ID	P02775
Sensitivity	2 pg/ml
Dilution Range	Dilution ratio of 1:500, concentration in serum and plasma is 50-150 ng/ml.
Expiration Date	Please enquire.



Human
CXCL7
PicoKine
ELISA Kit
standard ...