


Product Datasheet

Human NOV/CCN3 ELISA Kit (orb178689)

Description	Human NOV/CCN3 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: T32-M357
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
Clonality	Polyclonal
MW	39162 MW
Uniprot ID	P48745
Sensitivity	2 pg/ml
Expiration Date	Please enquire.



Human
NOV/CCN3
PicoKine
ELISA Kit
standa...