


Product Datasheet

Human CD13/Aminopeptidase N ELISA Kit (orb178692)

Description	Human CD13/Aminopeptidase N ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: K69-K967
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	109540 MW
Uniprot ID	P15144
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
CD13/Aminopeptidase
N PicoKine ELI...