


Product Datasheet

Human LAP(TGF-Beta1) ELISA Kit (orb178693)

Description	Human LAP(TGF-Beta1) ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L30-R278
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	44341 MW
Uniprot ID	P01137
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:100, concentration in serum and plasma is around 80 ng/ml.
Expiration Date	Please enquire.



Human
LAP(TGF-
beta1)
PicoKine
ELISA Kit ...