


Product Datasheet

Human soluble CD6 ELISA Kit (orb178696)

Description	Human soluble CD6 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: H18-E398
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	71801 MW
Uniprot ID	P30203
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
soluble CD6
PicoKine
ELISA Kit
sta...