


Product Datasheet

Mouse CD13/Aminopeptidase N ELISA Kit (orb178697)

Description	Mouse CD13/Aminopeptidase N ELISA Kit
Reactivity	Mouse
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: K69-S966
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	109651 MW
Uniprot ID	P97449
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:50, concentration in serum is around 200ng/ml
Expiration Date	Please enquire.



Mouse
CD13/Aminopeptidase
N PicoKine ELI...