


Product Datasheet

Human CCL23/MPIF-1 ELISA Kit (orb178699)

Description	Human CCL23/MPIF-1 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: R22-N120
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	13411 MW
Uniprot ID	P55773
Sensitivity	2 pg/ml
Expiration Date	Please enquire.



Human
CCL23/MPIF-
1 PicoKine
ELISA Kit
st...