


Product Datasheet

Mouse CCL9/MIP-1 gamma ELISA Kit (orb178700)

Description	Mouse CCL9/MIP-1 gamma ELISA Kit
Reactivity	Mouse
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: Q22-Q122
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	13871 MW
Uniprot ID	P51670
Sensitivity	1 pg/ml
Dilution Range	Dilution ratio of 1:1000, concentration in serum and plasma is 50-200 ng/ml.
Expiration Date	Please enquire.



Mouse CCL9
PicoKine
ELISA Kit
standard
C...