


Product Datasheet

Human TNFRSF14/HVEM ELISA Kit (orb178701)

Description	Human TNFRSF14/HVEM ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L39-V202
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	30392 MW
Uniprot ID	Q92956
Sensitivity	1 pg/ml
Dilution Range	Dilution ratio of 1:5, concentration in serum and plasma is around 800 pg/ml.
Expiration Date	Please enquire.



Human
TNFRSF14/HVEM
PicoKine ELISA
Kit s...