

Product Datasheet

Mouse TNFRSF14/HVEM ELISA Kit (orb178702)

Description	Mouse TNFRSF14/HVEM ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: Q29-V207
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	30171 MW
Uniprot ID	Q80WM9
Sensitivity	1 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 300-500 pg/ml.
Expiration Date	Please enquire.