


Product Datasheet

Band 4.1-like protein 1 EPB41L1 Antibody (orb178712)

Description

Band 4.1-like protein 1 EPB41L1 Antibody



WB analysis using Anti-Band 4.1-like pro...

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

A synthetic peptide corresponding to a sequence at the C-terminus of human EPB41L1 (585-604aa QDQERDTVFLKDNHLAIERK), different from the related mouse and rat sequences by one amino acid.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

98503 MW

Uniprot ID

[Q9H4G0](#)

Dilution Range

Western blot, 0.1-0.5µg/ml, Mouse, Rat, Human

Expiration Date

12 months from date of receipt.