

Product Datasheet

Anti-Mabthera monoclonal antibody (3D6H10) (orb1787827)

Description	Purified Mouse Monoclonal Antibody (Mab)
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Target	Mouse monoclonal antibody raised against Trastuzumab.
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	Mouse IgG
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.