

## Product Datasheet

### MS4A1/CD20 Antibody (C-term) (orb1787837)

|                            |                                                                                                                                                                               |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse Monoclonal Antibody (Mab)                                                                                                                                               |
| <b>Species/Host</b>        | Mouse                                                                                                                                                                         |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                                                             |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                  |
| <b>Tested Applications</b> | WB                                                                                                                                                                            |
| <b>Immunogen</b>           | 266-294 aa                                                                                                                                                                    |
| <b>Target</b>              | This MS4A1/CD20 antibody is generated from mice immunized with a KLH conjugated synthetic peptide between 266-294 amino acids from the C-terminal region of human MS4A1/CD20. |
| <b>Note</b>                | For research use only                                                                                                                                                         |
| <b>Isotype</b>             | IgG1                                                                                                                                                                          |
| <b>Clonality</b>           | Monoclonal                                                                                                                                                                    |
| <b>MW</b>                  | 33077 Da                                                                                                                                                                      |
| <b>Uniprot ID</b>          | <b>P11836</b>                                                                                                                                                                 |
| <b>NCBI</b>                | <b>NP_068769.2, NP_690605.1</b>                                                                                                                                               |
| <b>Dilution Range</b>      | WB: 1:3000                                                                                                                                                                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                               |