

Product Datasheet

HA Tag Antibody (orb1787856)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Target	KLH conjugated synthetic peptide encoding HA tag was used as antigen.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	45-50
Dilution Range	WB: 1:1000, WB: 1:1000, WB: 1:1000
Expiration Date	12 months from date of receipt.