

Product Datasheet

Erk1/2 Antibody (orb1787899)

Description	Mouse Monoclonal Antibody (Mab)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	28-240 aa
Target	Purified His-tagged Erk1/2 protein was used to produced this monoclonal antibody.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	43136 Da
Uniprot ID	P27361
Dilution Range	WB: 1:4000, FC: 1:25
Expiration Date	12 months from date of receipt.