

Product Datasheet

SNX1 Antibody (C-term) (orb1787935)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	426-455 aa
Target	This SNX1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 426-455 amino acids from the C-terminal region of human SNX1.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	59070 Da
Uniprot ID	Q13596
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.