

Product Datasheet

DPYSL5 Antibody (C-term) (orb1787987)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	478-512 aa
Target	This DPYSL5 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 478-512 amino acids from the C-terminal region of human DPYSL5.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	61421 Da
Uniprot ID	Q9BPU6
Dilution Range	IF: 1:25, WB: 1:2000
Expiration Date	12 months from date of receipt.