

Product Datasheet

Histone H3 Antibody (orb1788202)

Description	Purified Mouse Monoclonal Antibody (Mab)
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	1-156 aa
Target	This Histone H3 antibody is generated from a mouse immunized with Histone H3 recombinant protein.
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	15404 Da
Uniprot ID	P68431
Dilution Range	WB: 1:1000, WB: 1:3000
Expiration Date	12 months from date of receipt.