

Product Datasheet

ESR2 Antibody(C-term) (orb1788210)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	451-480 aa
Target	This ESR2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 451-480 amino acids from the C-terminal region of human ESR2.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	59216 Da
Uniprot ID	Q92731
NCBI	NP_001258805.1 , NP_001035365.1 , NP_001258806.1 , NP_001428.1 , NP_001201832.1 , NP_001201831.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.