

Product Datasheet

Phospho-ERBB2(Y877) Antibody (orb1788220)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT, WB
Immunogen	860-883 aa
Target	This ERBB2 Antibody is generated from rabbits immunized with a KLH conjugated synthetic phosphopeptide corresponding to amino acid residues surrounding Y877 of human ERBB2.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	137910 Da
Uniprot ID	P04626
NCBI	NP_004439.2, NP_001005862.1
Dilution Range	WB: 1:1000, DB: 1:500
Expiration Date	12 months from date of receipt.