

Product Datasheet

MAPK3/1 Antibody (Center) (orb1788245)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	195-228 aa
Target	This MAPK3/1 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 195-228 amino acids from the Central region of human MAPK3/1.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	43136 Da
Uniprot ID	P27361
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.