

Product Datasheet

UCHL1 Antibody (orb1788286)

Description	Mouse Monoclonal Antibody (Mab)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	1-243 aa
Target	This UCHL1 monoclonal antibody is generated from mouse immunized with UCHL1 recombinant protein.
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	24824 Da
Uniprot ID	P09936
NCBI	NP_004172.2
Dilution Range	IF: 1:25, WB: 1:1000, IHC-P: 1:25
Expiration Date	12 months from date of receipt.