

Product Datasheet

Phospho-HIST1H3B3(S10) Antibody (orb1788339)

Description

Affinity Purified Rabbit Polyclonal

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	8-35 aa
Target	This HIST1H3B3 Antibody is generated from rabbits immunized with a KLH conjugated synthetic phosphopeptide corresponding to amino acid residues surrounding S10 of human HIST1H3B3.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	15404 Da
Uniprot ID	P68431
NCBI	NP_003522.1, NP_003523.1, NP_066298.1, NP_003525.1, NP_003520.1, NP_003528.1, NP_003524.1, NP_003521.2, NP_003526.1, NP_003527.1
Dilution Range	WB: 1:500, WB: 1:500
Expiration Date	12 months from date of receipt.