

Product Datasheet

MDH2 Antibody (Center) (orb1788342)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	86-115 aa
Target	This MDH2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 86-115 amino acids from the Central region of human MDH2.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	35503 Da
Uniprot ID	P40926
NCBI	NP_001269333.1 , NP_001269332.1 , NP_005909.2
Dilution Range	WB: 1:1000, IHC-P: 1:10~50, FC: 1:10~50
Expiration Date	12 months from date of receipt.