

Product Datasheet

ADRA1D Antibody (N-term) (orb1788361)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	1-30 aa
Target	This ADRA1D antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 1-30amino acids from the N-terminal region of human ADRA1D.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	60463 Da
Uniprot ID	P25100
Dilution Range	WB: 1:1000, IHC-P: 1:25, IHC-P: 1:25, IHC-P: 1:25, FC: 1:25
Expiration Date	12 months from date of receipt.