

Product Datasheet

RAD51 Antibody (orb1788415)

Description	Purified Mouse Monoclonal Antibody (Mab)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Target	This RAD51 antibody is generated from a mouse immunized with a recombination protein from human.
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1, κ
Clonality	Monoclonal
MW	36966 Da
Uniprot ID	Q06609
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.