

Product Datasheet

CTSD Antibody (orb1788477)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Target	This antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between amino acids from human.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	44552 Da
Uniprot ID	P07339
Dilution Range	IF: 1:25, WB: 1:1000, IHC-P: 1:25, IHC-P: 1:25
Expiration Date	12 months from date of receipt.