

## Product Datasheet

### Functional BAFF-R (mouse) Antibody, mAb (orb1788507)

|                            |                                            |
|----------------------------|--------------------------------------------|
| <b>Description</b>         | Functional BAFF-R (mouse)<br>Antibody, mAb |
| <b>Reactivity</b>          | Mouse                                      |
| <b>Conjugation</b>         | Unconjugated                               |
| <b>Tested Applications</b> | FC                                         |
| <b>Concentration</b>       | 1mg/ml                                     |
| <b>Note</b>                | For research use only                      |
| <b>Clonality</b>           | Monoclonal                                 |
| <b>Uniprot ID</b>          | <b>Q9D8D0</b>                              |
| <b>Expiration Date</b>     | 12 months from date of receipt.            |