

## Product Datasheet

### p57KIP2 Antibody (orb1788897)

|                            |                                 |
|----------------------------|---------------------------------|
| <b>Description</b>         | p57KIP2 Antibody                |
| <b>Species/Host</b>        | Mouse                           |
| <b>Conjugation</b>         | Unconjugated                    |
| <b>Tested Applications</b> | IHC                             |
| <b>Note</b>                | For research use only           |
| <b>Clonality</b>           | Monoclonal                      |
| <b>Expiration Date</b>     | 12 months from date of receipt. |