


Product Datasheet

CD276 antibody (orb1791)

Description	Goat polyclonal antibody to CD276
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	B7-H3 / CD276
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 55kDa band observed in Human Placenta and Testis lysates (calculated MW of 57.2kDa according to NP_001019907.1). Recommended concentration: 0.1-3µg/ml.
Protein Sequence	KHSDSKEDDGQE
Clonality	Polyclonal
MW	57.2; 33.8
Entrez	80381
Dilution Range	ELISA: 1:8000, WB: 0.1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Placenta ...