

## Product Datasheet

### SLC6A4 Immunizing peptide (orb17948)

|                          |                                                                                               |
|--------------------------|-----------------------------------------------------------------------------------------------|
| <b>Description</b>       | SLC6A4 Immunizing peptide                                                                     |
| <b>Preservatives</b>     | 100µg of Lyophilised peptide, formulation<br>De-ionized (DI) water.                           |
| <b>Storage</b>           | Shipped at ambient temperature, store at -<br>20°C                                            |
| <b>Note</b>              | For research use only                                                                         |
| <b>Application notes</b> | This peptide was used as immunogen to<br>produce orb17948 (Serotonin Transporter<br>antibody) |
| <b>Protein Sequence</b>  | LHQGERETWGK                                                                                   |
| <b>Expiration Date</b>   | 12 months from date of receipt.                                                               |