



Product Datasheet

AGPS antibody (orb179854)


Description	Mouse monoclonal antibody to AGPS.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	recombinant human AGPS (amino acids 158-384)
Target	AGPS
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml. Intracellular staining.
Isotype	Mouse IgG2a
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	O00116
Entrez	8540
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml. Intracellular staining. Western blotting: Recommended dilution: 1-2 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of
K562 cells
(red-filled)
fr...



Flow
cytometry
intracellular
staining pa...



Western blotting analysis of human AGPS